

RPT0035LQ

Leader in Biomolecular Solutions for Life Science



Recombinant His-TrxA-SUMO Protein

Catalog No.: RPT0035LQ

Recombinant

Sequence Information

Species	Gene ID	Swiss Prot
<i>E. coli</i>	/	/

Tags

/

Synonyms

Product Information

Source	Purification
	≥ 95 % as determined by SDS-PAGE.

Endotoxin

Please contact us for more information.

Formulation

0.22 µm filtered solution of PBS, pH 7.4.

Reconstitution

Background

Basic Information

Description

Recombinant His-TrxA-SUMO Protein is produced by *E. coli* expression system.

Bio-Activity

Shipping

Lyophilized proteins are shipped at ambient temperature.
Liquid proteins are shipped at frozen temperature with dry ice.
Upon receipt, store it immediately at the temperature recommended below.

Operational Notes

For your safety and health, please wear a lab coat and disposable gloves for handling.

Storage

Store at -80°C
Avoid repeated freeze/thaw cycles.


Contact



www.abclonal.com

* For your safety and health, please wear a lab coat and disposable gloves when handling.

Validation Data



Recombinant His-TrxA-SUMO
Protein was determined by SDS-
PAGE under reducing conditions with
Coomassie Blue.