

YR0117

Leader in Biomolecular Solutions for Life Science



Human CD33 Monoclonal Antibody, Endotoxin 0.05 EU/mg

Catalog No.: YR0117

Basic Information

Molecular Weight

150 kDa

Endotoxin

<0.05EU/mg
(<0.00005EU/μg)Determined by LAL
gel clotting assay

Sterility

0.2 μm filtration

Aggregation

<5% Determined by SECP

Purity

>95% Determined by SDS-PAGE

Background

Enables protein phosphatase binding activity and sialic acid binding activity. Involved in several processes, including negative regulation of cytokine production; negative regulation of monocyte activation; and positive regulation of protein tyrosine phosphatase activity. Located in several cellular components, including Golgi apparatus; external side of plasma membrane; and peroxisome.

Reported Applications

ELISA,neutralization,functional assays such as bioanalytical PK and ADA assays,and those assays for studying biological pathways

Immunogen Information

Clone

Gemtuzumab Biosimilar

Isotype

Human IgG4 kappa

Immunogen

Human CD80

Recommended Isotype Control(s)

In Vivo Grade Recombinant Human IgG4-S228P Kappa Isotype Control Antibody

Recommended Dilution Buffer

1×PBS pH 7.3

Contact



www.abclonal.com

Product Information

Production

Purified from cell culture supernatant in an animal-free facility

Purification

Protein A or G purification

Storage

Store at -20°C