

YR0118

Leader in Biomolecular Solutions for Life Science



Human CTLA-4 Monoclonal Antibody

Catalog No.: YR0118

Basic Information

Molecular Weight

150 kDa

Endotoxin

<1EU/mg (<0.001EU/μg) Determined by LAL gel clotting assay

Sterility

0.2 μm filtration

Aggregation

<5% Determined by SECP

Purity

>95% Determined by SDS-PAGE

Background

Ipilimumab is a human cytotoxic T-lymphocyte antigen 4 (CTLA-4) blocking antibody used to treat metastatic or unresectable melanoma. Ipilimumab is a fully humanized IgG1 monoclonal antibody that blocks cytotoxic T lymphocyte antigen-4 (CTLA-4). Blocking CTLA-4 removes an inhibitory signal from reducing the activity of T lymphocytes. Cytotoxic T-lymphocyte antigen-4 (CTLA-4) is an inhibitory molecule that competes with the stimulatory CD28 for binding to B7 on antigen presenting cells. CTLA-4 and CD28 are both presented on the surface of T-cells. Ipilimumab is a human IgG1 that binds CTLA-4, preventing the inhibition of T-cell mediated immune responses to tumors.

Reported Applications

ELISA, neutralization, functional assays such as bioanalytical PK and ADA assays, and those assays for studying biological pathways

Immunogen Information

Clone

Ipilimumab Biosimilar

Isotype

Human IgG1 kappa

Immunogen

Human CTLA-4

Recommended Isotype Control(s)

In Vivo Grade Recombinant Human IgG1 Kappa Isotype Control Antibody

Recommended Dilution Buffer

1×PBS pH 7.0

Contact



www.abclonal.com

Product Information

Production

Purified from cell culture supernatant in an animal-free facility

Purification

Protein A or G purification

Storage

Store at 2 - 8°C. 2 - 8°C for up to 4 weeks and -80°C for long term storage (Avoid repeated freezing and thawing)