

YR0503

Leader in Biomolecular Solutions for Life Science



# Human CD123 Monoclonal Antibody ( N297Q)

Catalog No.: YR0503

## Basic Information

---

### Molecular Weight

150 kDa

### Endotoxin

<1EU/mg (<0.001EU/μg) Determined by LAL gel clotting assay

### Sterility

0.2 μm filtration

### Aggregation

<5% Determined by SECP

### Purity

>95% Determined by SDS-PAGE

## Background

---

Pivekimab is a human IgG1 monoclonal antibody that targets interleukin-3 (IL-3) and CD123. Pivekimab is used to synthesize pivekimab sunirine, a CD123-directed antibody-active molecule conjugate (ADC). Pivekimab is used in the study of blastic plasmacytoid dendritic cell neoplasm (BPDCN).

## Reported Applications

---

ELISA, neutralization, functional assays such as bioanalytical PK and ADA assays, and those assays for studying biological pathways

## Immunogen Information

---

### Clone

Pivekimab Biosimilar

### Isotype

Human IgG1 kappa

### Immunogen

### Recommended Isotype Control(s)

### Recommended Dilution Buffer

1×PBS pH 7.3

## Contact

---



[www.abclonal.com](http://www.abclonal.com)

## Product Information

---

### Production

Purified from cell culture supernatant in an animal-free facility


### Purification

Protein A or G purification


### Storage

Store at 2 - 8°C. 2 - 8°C for up to 4 weeks and -80°C for long term storage (Avoid repeated freezing and thawing)

## Validation Data



Direct ELISA binding curve demonstrating the recognition of Human Anti-Human CD123 (Research Grade Pivekimab (N297Q) Biosimilar) Monoclonal Antibody to CD123. The target protein was coated onto the microplate well surface, followed by binding of the antibody. A donkey anti-human IgG HRP conjugate was used for detection.



Human CD123 SPR assay. Determined through SPR assay, the Human Anti-Human CD123 (Research Grade Pivekimab (N297Q) Biosimilar) Monoclonal Antibody is Capable of binding to Human Trop-2 with an affinity constant of 29.26pM.